

Verifying Connection to Spot

Verifying Connection:

1. Make sure Tape Measure and your device are on the same network.
2. Ping Spot by opening command prompt and paste:

Windows 10/11:


```
python.exe -m bosdyn.client SPOT-IP id # Replace SPOT-IP with the network address given to it
```

For Mac & Ubuntu users:

```
python3 -m bosdyn.client SPOT-IP id # Replace SPOT-IP with the network address given to it
```

This should return Tape Measure's name and ID.

*Should look something like this:



```
Command Prompt
Microsoft Windows [Version 10.0.22621.2283]
(c) Microsoft Corporation. All rights reserved.

C:\Users\caich>python.exe -m bosdyn.client 192.168.200.39 id
spot-BD-31100006      000060208679A4000029 K9: Tape Measure spot (V3)
Software: 3.3.0 (fe5d7bf59da75a371c238e6d7fbc321e84918e40 2023-05-28 11:50:08)
Installed: 2023-05-28 12:55:42

C:\Users\caich>
```

Revision #1

Created 22 January 2024 03:36:06 by Dwight Howard, II

Updated 22 January 2024 03:47:39 by Dwight Howard, II